

PRODUCT	DESCRIPTION AND USAGE	PRICE IN EURO [€]
SIGMA-1SA	<p>LEYCOM® SIGMA-1SA: (single segment) For the measurement and display of intra-ventricular volume and pressure in mice. (Compatible with Millar and Scisense Mouse Catheters)</p> <p><u>Included are the following accessories:</u></p>	9,950.00
CFL-CONNT-CD	<p>CD LEYCOM Conduct NT (M) software, version: V....</p> <p><i>Usage: Windows based Software Package, for the acquisition, display, analysis and storage of pressure and volume data</i></p>	
1SA-CONNT-MAN	LEYCOM® Conduct NT manual , for Sigma-1SA, version: V....	44.00
CCBL-BRM-0.5M	<p>Catheter Connection Cable, Binder Redell Mouse, Single Field, length 0.5 meter</p> <p><i>Usage: to connect a Millar mouse catheter to the Sigma-1SA together with MILLAR-TEC10D (to be supplied by Millar)</i></p>	216.00
1SA-USB-1.0M	<p>USB Cable, length 1.0 meter</p> <p><i>Usage: Cable to connect the USB-port of the computer to the Sigma-1SA</i></p>	27.50
SIPWRCRD-XXX	<p>Powercord for Leycom® Sigma-1SA <i>* Country specific, see below</i></p>	73.50
* SIPWRCRD-EUR / SIPWRCRD-USA / SIPWRCRD-UK /	* Powercord Europe for Leycom® Sigma-1SA Powercord United States for Leycom® Sigma-1SA / Powercord United Kingdom for Leycom® Sigma-1SA DF /	

This price list supercedes all other price lists before January 18, 2007.

Prices are valid till December 31 2007 for those countries in the Eurozone which do not have local representation.

Prices are indicative for countries outside the Eurozone and which do not have local representation.