

PRODUCT	DESCRIPTION AND USAGE
	4 French catheters without lumen:
CA-41063-PN	CD Leycom Pressure/Volume Combination Catheter: 4 French, 12 electrodes, 1 pressure sensor between electrodes 5 and 6, with pigtail, no lumen <i>Spacing: 6 mm. spacing between the electrodes</i>
CA-41103-PN	CD Leycom Pressure/Volume Combination Catheter: 4 French, 12 electrodes, 1 pressure sensor between electrodes 5 and 6, with pigtail, no lumen <i>Spacing: 10 mm. spacing between the electrodes</i>
	7 French catheters without lumen:
CA-71083-PN	CD Leycom Pressure/Volume Combination Catheter: 7 French, 12 electrodes, 1 pressure sensor between electrodes 5 and 6, with pigtail, no lumen <i>Spacing: 8 mm. spacing between the electrodes</i>
CA-71103-PN	CD Leycom Pressure/Volume Combination Catheter: 7 French, 12 electrodes, 1 pressure sensor between electrodes 5 and 6, with pigtail, no lumen <i>Spacing: 10 mm. spacing between the electrodes</i>
CA-71123-PN	CD Leycom Pressure/Volume Combination Catheter: 7 French, 12 electrodes, 1 pressure sensor between electrodes 5 and 6, with pigtail, no lumen <i>Spacing: 12 mm. spacing between the electrodes</i>
	7 French catheters with lumen:
CA-71083-PL	CD Leycom Pressure/Volume Combination Catheter: 7 French, 12 electrodes, 1 pressure sensor between electrodes 5 and 6, with pigtail, with 0.025" guide wire lumen <i>Spacing: 8 mm. spacing between the electrodes</i>
CA-71103-PL	CD Leycom Pressure/Volume Combination Catheter: 7 French, 12 electrodes, 1 pressure sensor between electrodes 5 and 6, with pigtail, with 0.025" guide wire lumen <i>Spacing: 10 mm. spacing between the electrodes</i>
CA-71123-PL	CD Leycom Pressure/Volume Combination Catheter: 7 French, 12 electrodes, 1 pressure sensor between electrodes 5 and 6, with pigtail, with 0.025" guide wire lumen <i>Spacing: 12 mm. spacing between the electrodes</i>
	7 French catheters without lumen, with double pressure-sensor:
CA-72123-PN	CD Leycom Pressure/Volume Combination Catheter: 7 French, 12 electrodes, 1 pressure sensor between electrodes 5 and 6, 1 pressure sensor 50 cm. proximal from electrode 12, with pigtail, no lumen <i>Spacing: 12 mm. spacing between the electrodes</i>
CA-72103-PNA	CD Leycom Pressure/Volume Combination Catheter: 7 French, 12 electrodes, 1 pressure sensor between electrodes 5 and 6, 1 pressure sensor 2 cm. proximal from electrode 12 (in the aorta), with pigtail, no lumen <i>Spacing: 10 mm. spacing between the electrodes</i>
	6 French pressure-only catheters: (no FDA-approval)
CA-61000-PL	CD Leycom Pressure Only Catheter: 6 French, 1 pressure sensor at distal tip, with pigtail, 0.028" guide wire lumen, with side holes, max. pressure 1000 P.S.I., max. bolus 35 ml., flow rate 15 ml./sec., straight
CA-61000-PLB	CD Leycom Pressure Only Catheter: 6 French, 1 pressure sensor at distal tip, with pigtail, 0.028" guide wire lumen, with side holes, max. pressure 1000 P.S.I., max. bolus 35 ml., flow rate 15 ml./sec., 155° angle 6 centimeter from distal tip